

封面

Operation

1. Download App

To download mobile application to your Mobile device:

• Android: Go to Google play download "OPCOM Care" App.

• iOS: Go to App store download "OPCOM Care" App.

2. Paring Device to Mobile Device :

OPCOM Care unit to Mobile Device:

• Turn on Wi-Fi on your mobile device, connect to Duck_xxxxxxxxxx from the Wi-Fi list.

(Option)USB On To Go (OTG) Mode :

• Plug Type c into Type c port at the bottom of OPCOM Care Unit.

• Plug another head of Type c into OPCOM Care System Tablet.

Care

USB OTG

10.6" Tablet

3. Live View :

Open application from your Mobile or OPCOM Care system tablet.

4. Remark :

Mobile / Application User's Manual available at: www.opcomcare.com

Weatherproof

USB OTG

OS + APP

Windows PC AP

USB Charger

6mm Side View

Microview Scope122°

Disposable Cover

For detailed warranty information,Please visit www.opcomcare.com

OPCOM

Duck

Intra-Oral

Camera

Quick Guide

CL481D

OPCOM Duck

A personal health scope to check yours and your loved ones ears, nose,eyes, skin and, hair.With its built-in mini 960P camera in the 6mm head, pictures and videos can be captured for real-time viewing, or video that can be stored in your personal files and forwarded to a medical professional for review and diagnosis. Easy to use on small children and infants or pets....

1

2

3

4

5

6

7

Front View

Side View

8

(Option) Link To Tablet

1. EndoCam Camera

960P camera.

2. Snapshot button

• Press button to capture image.

• Press button for 1 second to take a snapshot picture.

• Hold button down for 2-3 seconds to start video recording.

3. Brightness button

Quickly press button and release it to increase / decrease light.

4. Light indicator

No connect to USB cable

1. Power on: ● blink 10s link to WiFi mode, then keep ● on

2. Power off: ● blink >3s to exit WIFI mode

3. Low power: ● continue blinking

4. Photograph: ● blink once

5. Recording: ● continue continually

Stop recording: keep ● on

6. EndoCam camera cconnection not good: ●● alternate blinking, then power off

Connect to USB cable:

1. Charging mode: ● on

2. Fully charged: ● on

5. Power ON/OFF Button

PC Mode:

Press <3 second to power ON and PC Mode

Press 1 second to power OFF.

WIFI Mode:

Press >3 second to WIFI ON.

Press 3 second to WIFI OFF.

6. EndoCam Camera Head

Camera and LED light module.

7. Battery

Location where rechargeable battery is placed into the unit .

8. USB Cable

Type C port at the bottom of the unit.

9. Charge

Support operating time up to 7.5 hours per each full charge .

FCC STATEMENT

1. This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:
(1) This device may not cause harmful interference, and
(2) This device must accept any interference received, including interference that may cause undesired operation.
2. Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.
NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
Reorient or relocate the receiving antenna.
Increase the separation between the equipment and receiver.
Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
Consult the dealer or an experienced radio/TV technician for help.

尺寸 : 140*140mm 材质 : 100g铜板

P414811401G1-481D说明书-20190408

C 50

Y 100